

PE CARE PORTABLE ALUMINIUM SCOOTER WHEELCHAIR RAMPS - SWL 430KG



SKU: PEC-8330

SPECIFICATIONS

Closed Length - 74.5cm

Extended Length - 150cm

Width - 20cm

Material - Aluminum

Weight per ramp: 6kg each

Total Weight: 12kgs

WARRANTY

1 Year Structural Warranty

PRODUCT DESCRIPTION

The PE Care Portable Aluminium Scooter Wheelchair Ramps are the perfect accessory for your Mobility Scooter, Wheelchair or Power Chair. They make transporting your Scooter much easier enabling you to drive your scooter straight up into a vehicle or onto a trailer. Being portable, they are also handy for around the home. You can quickly set them up to get your mobility aid in and out of your front or back door, or up steps.

The ramps have a telescopic design. You can adjust their length from 74.5cm to 150cm with a width of 20cm. The ramps have a non slip top surface to help the wheels of your mobility aid grip to the ramp.

These Ramps are sold in pairs and have handles to make them easier to carry. This makes for easy transportation and storage. Being just under 75cm, they can fit easily in the boot of your car; or you can store them standing up behind a door or underneath the bed.