

## [SHOPRIDER 889D EXTRA WIDE SEAT MOBILITY SCOOTER](#)



**SKU:** SHO-8828

### **SPECIFICATIONS**

Overall Length: 148 cm

Overall Width: 91 cm

Overall Height: 98 cm

Turning Circle: 255 cm

Speed: 10km/h

Maximum range: Up to 40km

Safe climb angle: 9°

Maximum climb angle: 12°

Motor: 2.0Hp high torque motor (2 Pole)

Batteries: 2 x 75 Amp minimum (supplied)

Charger: 8 Amp



Tyres (Pneumatic): Ø330mm – Front + Rear

Product Weight: 159kg – with batteries

Maximum User Weight: 225kg – on flat ground

Please note: Single person use only

## **WARRANTY**

### **Shoprider Warranty**

Terms and Conditions

#### **Structural Frame - 10 Year Warranty:**

Shoprider Australia warrants that the frame will be free from defects under normal use, for a period of ten (10) years from the date of purchase by the original customer. Shoprider Australia will repair or replace, at our option, free of charge any part found upon examination by Shoprider Australia to be defective in material and/or workmanship within the warranty period. All parts claimed under warranty must be returned to Shoprider Australia in their present condition via the original place of purchase (SHOPRIDER AUTHORISED DEALER). Only an authorised "WARRANTY DEALER" of Shoprider Australia, following approval by Shoprider Australia, can undertake warranty work.

#### **Drive Train Assembly - 18 Month Warranty:**

Shoprider Australia warrants that the drive train assembly will be free from defects under normal use, for a period of eighteen (18) months from the date of purchase by the original customer. Shoprider Australia will repair or replace, at our option, free of charge any part found upon examination by Shoprider Australia to be defective in material and/or workmanship within the warranty period. All parts claimed under warranty must be returned to Shoprider Australia in their present condition via the original place of purchase (SHOPRIDER AUTHORISED DEALER). Only an authorised "WARRANTY DEALER" of Shoprider Australia, following approval by Shoprider Australia, can undertake warranty work. Drive Train Components Include:

Transaxle, Motor and Magnetic Brake.

#### **OTHER PARTS - 12 Month Warranty:**

Shoprider Australia warrants other parts will be free from defects under normal use, for a period of twelve (12) months from the date of purchase by the original customer. Shoprider Australia will repair or replace, at our option, free of charge **any part** found upon examination by Shoprider Australia to be defective in material and/or workmanship within the warranty period. All parts claimed under warranty must be returned to Shoprider Australia in their present condition via the original place of purchase (SHOPRIDER AUTHORISED DEALER). Only an authorised "WARRANTY DEALER" of Shoprider Australia, following approval by Shoprider

Australia, can undertake warranty work.**Exclusions:**

1. All transportation costs associated with the submission and return of parts for repair or replacement under warranty (as determined by Shoprider Australia) are the responsibility of the original customer.
2. Normal wear items are not covered under this warranty. These include such items as brushes, tyres, tubes, brake pads, upholstery, seating and arm rests.
3. Damage caused by:
  - >Excessive, abusive, misuse, negligence or accident.
  - >Improper maintenance, maintenance carried out by non-approved "WARRANTY DEALER".
  - >Improper operation or storage.
  - >Any Commercial use.
  - >Fitment of NON-Manufacturer approved accessories.Resulting damage attributed to the NON-Manufacturer approved accessory or correction of the unit to operate as intended once the accessory have been removed will not be covered under warranty.
4. Where an authorised dealer has not lodged all the required documentation at the appropriate time, the warranty is void.
5. All warranties provided do not extend to purchasers of used or second hand equipment.

## **PRODUCT DESCRIPTION**

The Shoprider 889D Extra Wide Seat Mobility Scooter has been specifically designed to assist with the growing demand for an extra-wide 4 wheel mobility scooter. Its seat spans over 80cm with a max user weight of 225kg on flat ground. The body of the 889D Scooter has an overall width of 91cm, accommodating for this later seat design and still keeping the scooter stable and balanced to the ground. Its adjustable front and rear suspension provides added comfort and support for all types of users. The large 2.0Hp high torque rear wheel-driven motor along with two 75amp batteries gives it a generous 40km driving range and a top speed of 10km/h. The 889D Gopher comes complete with 13" pneumatic tyres, perfect for traveling over various types of terrain. Available in red only and comes with a 3-year structural warranty.

Additional features of the Shoprider 889D Scooter are:

- Gas Adjustable Steering Tiller
- Loop Handles With Easy Finger Controls
- High / Low Speed Settings
- Electromagnetic Brakes
- Twin Headlights

- Rear Lights
- Indicators
- Hazard Lights
- Forward And Backwards Sliding Seat
- Front Shopping Basket
- Easy Lift Armrests.

Please note: Single person use only

**Please note: We do not ship to all areas in NSW. For more information, please contact one of our friendly staff on 1300 108 280 or [sales@mobilityandwellness.com.au](mailto:sales@mobilityandwellness.com.au)**