

SHEAR COMFORT NEONATAL AND PAEDIATRIC OVERLAY



SKU: SHE-6699

SPECIFICATIONS

Size: Neonatal: 50 x 30cm

Paediatric: 74 x 45cm

Colour: White And Blue

Material: XD1900 Fabric-Backed Medical Sheepskin

Code: Neonatal - White - XD1900 - 005221

Neonatal - Blue - XD1900 - 002251

Paediatric - White - XD1900 - 003135

Washing Instructions

Machine washable up to 80 degrees Celsius

The fabric-backed material on which the wool is woven makes it easy to wash up to 80 degrees Celsius - for a hygienic wash.

Shear Comfort XD1900 does not require special shampoo, and can be washed with normal household detergent up to a temperature of 80°C, although you can use Shear Comfort Wool Care Shampoo to condition the wool and utilise the fungicidal and mitricidal properties of the shampoo.

50 washes guaranteed

Shear Comfort XD1900 products have a warranty of up to 50 washes when following the manufacturer's washing instructions.

Drying Instructions

To air-dry Shear Comfort footwear or targeted protectors, release all fastenings and place fleece-side down on a drying rack and away from direct sunlight or heat sources.

Do not iron

Do not bleach

WARRANTY

One year warranty from date of purchase, T&Cs apply

VARIATIONS



Image	SKU	Size	Colour
	SHE-6699-9830	Neonatal Overlay	White

Image	SKU	Size	Colour
	SHE-6699-9831	Paediatric Overlay	White

PRODUCT DESCRIPTION

The Shear Comfort Neonatal & Paediatric Overlay provides comfort for babies and toddlers delicate skin.

The Neonatal Overlay is the perfect size for newborns as it is small enough size-wise to place in an incubator, cot or Moses basket to protect your baby's skin and help them sleep soundly.

For older babies or toddlers, the Paediatric Overlay fits perfectly in a pushchair, cot bed or even on the floor for a cosy patch of Sheep's wool to rest and play on.

There have been a number of studies conducted that illustrate the benefits for babies sleeping on and under wool. The rate of weight gain in underweight newborns was 61% higher when sleeping on a wool overlay compared with a cotton sheet. Jaundiced newborns sleeping on wool were more settled when sleeping compared with those on cotton. Only around 30% of babies on wool cried compared with 67% on cotton.

Made from XD1900, this is medical sheep's wool unique to Shear Comfort, with a fibre density higher than even the very best natural sheepskins - which gives it superior pressure redistribution properties. Providing the ideal microclimate for a good sleep, looking after temperature and perspiration regulation, with the excellent air redistribution properties of a wool product.

Made with natural properties of fabric-backed medical sheepskin wool, unique to Shear Comfort, ensures that comfort and protection are provided. The moisture dissipation qualities of the wool ensure any perspiration is effectively wicked away from the skin.

All Shear Comfort products are machine washable and are guaranteed against deterioration for up to 50 washes. Both styles are available to purchase in white and blue. The neonatal measures 50 x 30cm and the Paediatric measures 74 x 45cm.

Please note: As these are hand made sometimes there is a small dispatch delay of 1-2 weeks as they need to be made up.

